



**Supermed**

Globally, Sustainable

# Personal Protective

product catalog



## ***All Gowns Series***

1. Surgical Gown
2. Isolation Gown
3. CPE Gown
4. Disposable PE Apron



## ***Scrub Suits & Coats Series***

1. PP/SMS Scrub Suits
2. SMS Jacket
3. PP/SMS Lab Coat



## ***Coveralls Series***

1. TYPE 3B/4B/5B/6B Protective Coverall
2. TYPE 4B/5B/6B SF Coverall With Reflective Stripes
3. TYPE 4B/5B/6B SF Coverall With Tape
4. TYPE 4/5/6 SF Coverall & TYPE 5/6 SMS Coverall
5. TYPE 5/6 SMS Coverall FIRE RETARDANT
6. PP/SMS/PP+PE Coverall With Hood



## ***Face Masks Series***

1. Surgical Face Mask
2. Disposable Children Face Mask
3. FFP2 Folded Respirator
4. FFP2 Duckbill Respirator
5. FFP3 Valved Respirator
6. KN95 Respirator
7. KF94 Face Mask
8. Disposable Dust Mask
9. Disposable Face Shield



## **Caps Series**

1. Bouffant/Round/Nurse Cap
2. Disposable Surgical/Doctor Cap
3. Disposable Surgical/Doctor Cap
4. Mob Cap/Clip Cap



## **Gloves Series**

1. Latex Gloves & Nitrile Gloves
2. Vinyl PVC Gloves & PE Gloves



## **Shoe Covers, Boots & Sleeves Series**

1. Disposable Non-skid PP Shoe Cover
2. Disposable PP Shoe Cover
3. Disposable PP+CPE Shoe Cover
4. Disposable PP/PE/PP+PE Boot Covers & Microporous Boot Cover With PVC Sole
5. Microporous Boot Cover With Tapes
6. Sleeves



## **Bed Sheets Series**

1. Disposable PP/PP+CPE/SMS Bed Cover
2. Disposable PP/PP+CPE/SMS Bed Sheets



*Easy Life*

# All Gowns

series



## Surgical Gown

- » 35~70gsm SMS/SMMS non-woven fabric
- » 68gsm BVB material
- » 60gsm spunlace non-woven fabric + breathable composite film
- » With long sleeves, knitted cuffs, velcro on the neck back, ties at the waist, paper card on the tie, with/without reinforcement, with/without wrap, with/without hand towels, EO sterilization or not
- » Size: 140\*130cm/150\*130cm/155\*130cm etc
- » Customizable according to customer needs





## Isolation Gown

- » PP: 15~50gsm
- » PP+PE: 18~70gsm
- » SMS/SMMS: 25~60gsm
- » Microporous: 45~65gsm
- » Long sleeves, ties/velcro on the neck, ties on waist, elastic/knitted cuff, sewing or ultrasonic
- » Color: white, yellow, blue, green, pink, purple etc



## CPE Gown

- » CPE film
- » Apron style, long sleeves, back empty, thumb loops/elasticated wrists, 2 ties
- » Color: blue, yellow etc



## Disposable PE Apron

- » Size: 28"x42", 28"x46", 28"x48", 28"x50", 28"x53"
- » Thickness: PE film 0.016~0.10mm





*Easy Life*

# Scrub Suits & Coats

series






## PP / SMS Scrub Suits

- » 20~55gsm
- » short sleeves, V neck or round neck, and different pockets are available on shirts/trousers
- » Color: dark blue, green, pink, purple, etc

## SMS Jacket

- » 25~60gsm
- » long sleeves, knitted collar and wrists, with/without pocket, 5 buttons
- » Color: dark green, dark blue, red, purple etc



## PP/SMS Lab Coat

- » 20~60gsm
- » long sleeves, shirt collar/knitted collar, with/without pocket, 5 buttons/zipper/velcro, open/elastic/knitted wrists
- » Color: white, green, blue, red, purple etc





*Easy Life*

# Coveralls

series



## Type 3B/4B/5B/6B Protective Coverall

- » 90gsm SMS+PE fabric
- » With Taped Seam, with 3 pieces hood and boots, double zipper system covered by adhesive flap, elastic at cuffs and waist; heat sealed seams by tapes
- » Color: white/blue/yellow etc
- » Size: S-5XL



## TYPE 4B/5B/6B SF Coverall With Reflective Stripes

- » 65gsm Microporous breathable fabric
- » With reflective stripes, with hood, with zipper covered by adhesive flap, elastic at cuffs, ankles and waist
- » Color: white/blue/yellow/red etc
- » Size: S-5XL





## **TYPE 4B/5B/6B SF Coverall With Tape**

- » 65gsm Microporous breathable fabric
- » With Taped Seam, with 3 pieces hood and boots, with zipper covered by adhesive flap, elastic at cuffs, ankles and waist; heat sealed seams by tapes
- » Color: white/blue/yellow/red etc
- » Size: S-5XL





## TYPE 4/5/6 SF Coverall

- » 55~60gsm Microporous breathable fabric
- » With hood, with zipper covered by adhesive flap, elastic cuffs, ankles and back waist. bound seams
- » Color: white/blue/yellow etc
- » Size: S-5XL

## Type5/6 SMS Coverall

- » 45~65gsm
- » SMS fabric with hood, with zipper, covered by adhesive flap, elastic at hood, cuffs, ankles and waist. serged seam inside
- » Color: white/blue/yellow etc
- » Size: S-5XL

## TYPE 5/6 SMS Coverall FIRE RETARDANT

- » 55gsm SMS fabric
- » With hood, with zipper, covered by adhesive flap, elastic cuffs, ankles and waist. serged seam inside
- » Color: white/blue/green/yellow etc
- » Size: S-5XL



## PP/SMS/PP+PE Coverall With Hood

- » 25~70gsm
- » With hood, elastic cuffs, ankles and waist, with zipper
- » Color: white, blue, yellow, dark blue etc
- » Size: S-5XL





*Easy Life*

# Face Masks

series



## Surgical Face Mask

- » PP/meltblown/PP material
- » 3-ply, BFE $\geq$ 99%, with ties or with earloop (round elastic/band), adjustable integrate nose bar, latex-free
- » Size: 17.5\*9.5cm
- » Color: white/blue/green/pink/yellow/black etc
- » Packing: 50pcs/box, 2000pcs/ctn



## Disposable Children Face Mask

- » 3-ply, with earloop, adjustable integrate nose bar, latex-free
- » Size: 14.5\*9.5cm
- » Color: white/blue/green/pink/yellow/black printed etc
- » Packing: 50pcs/box, 2000pcs/ctn
- » Customizable according to customer needs



## FFP2 Folded Respirator

- » FFP2 respirator
- » Flat-fold, with or without valve
- » Color: white



## FFP2 Duckbill Respirator

- » Duckbill shape
- » Color: white
- » Packing: 20pcs/box, 400pcs/ctn







## FFP3 Valved Respirator

- » Color: white
- » Packing: 20pcs/box, 200pcs/ctn

## KN95 Respirator

- » 5ply including double layers meltblown material, falt-fold, with elastic ear band, or with headworn elastic band. with adjustable nose bar
- » Color: white/blue/black etc
- » KN95 is available packed in bags or boxes

## KF94 Face Mask

- » 4-ply, with earloop
- » Size: 20.5\*8cm
- » Color: white/blue/green/pink/yellow/black etc





## Disposable Dust Mask

- » 120~220gsm PP material, or 200~280gsm terylene material
- » Moulded cup shape, with elastic earloop. with aluminum nose bar
- » Color: white/blue etc

## Disposable Face Shield

- » PET shield with sponge and elastic band
- » Size: 33x22cm
- » packing: 10pcs/bag, 200pcs/ctn





*Easy Life*

# Caps

s e r i e s





## Bouffant/Round/Nurse Cap

- » PP/SMS/Spunlace/Microporous
- » 10~40gsm fabric
- » Hand-made, sewn with elastic
- » Color: white/blue/green/re/black etc





## Disposable Surgical/Doctor Cap

- » 25~40gsm
- » Machine-made, elastic or ties on back
- » Color: blue/green etc





## Disposable Surgical/Doctor Cap

- » 20~40gsm
- » Hand-made, ties on back
- » Color: blue/green/white/diamond pattern etc
- » Soft absorbent sweat pad on forehead is available



## Mob Cap/Clip Cap

- » 10~16gsm PP fabric
- » Machine-made, single/double elastic
- » Color: white/blue/green/yellow/red/black etc



*Easy Life*

# Gloves

s e r i e s





## **Nitrile Gloves**

- » Size: S/M/L/XL
- » Color: white/blue/green/black etc

## **Latex Gloves**

- » Size: S/M/L/XL
- » Color: white



## Vinyl PVC Gloves

- » Size: S/M/L/XL
- » Color: white/green/black etc

## PE Gloves

- » Size: S/M/L/XL
- » Color: clear etc
- » Long PE gloves are available



*Easy Life*

# Shoe Covers, Boots & Sleeves

s e r i e s



## Disposable Non-Skid PP Shoe Cover

- » 30~40gsm PP material
- » Machine-made, elastic on top and bottom
- » Color: white/blue/green/pink etc





## Disposable PP Shoe Cover

- » 25~40gsm PP material
- » Machine-made, elastic on top and bottom
- » Color: white/blue/green/pink etc
- » With Conductive Strip is available



## Disposable PP+CPE Shoe Cover

- » Machine-made, elastic on top
- » Color: white PP + white/blue CPE
- » Size: 15x36cm, 15x41cm





## Disposable PP/PE/PP+PE Boot Covers

- » Elastic on top, with/without ties
- » Color: white/yellow/green etc
- » Size: 48x38cm, 30x38cm



## Microporous Boot Cover With PVC Sole

- » 45~65gsm microporous
- » elastic on top and ankles, with PVC sole, with/without ties
- » Color: white etc

## Microporous Boot Cover With Tapes

- » 45~65gsm microporous
- » heat sealed seams with tap, elastic on top, with tie straps
- » Color: white+ white/blue/red tape



## Sleeves

- » Elastic/knitted cuff
- » 20~40gsm PP
- » 40~65gsm Microporous
- » 40~65gsm fabric
- » PE film
- » Color: white/blue





*Easy Life*

# Bed Sheets

series





## Disposable PP/PP+CPE/SMS Bed Cover

- » Elastic on 4 corners, or all round the edges
- » Size: 80\*180cm, 90\*200cm
- » Color: white/blue/pink etc





## Disposable PP/PP+CPE/SMS Bed Sheets

- » Size: 80\*180cm, 90\*200cm
- » Color: white/blue/pink etc



**Super<sup>+</sup>med**

Globally, Sustainable

 **+65 9055 9545**

 **info@supermed.com.sg**

 **www.supermed.com.sg**